

litilit

BIOLIT 2

Femtosecond Fiber Laser for Biophotonics
1050 nm, 70 fs, 2 W, 15 – 40 MHz



WWW.LITILIT.COM

CLEAN PULSES LEAD TO SHARP IMAGES

FEATURES

- Ultra-short and clean pulses
- Robust and stable
- Flexible repetition rate
- Maintenance-free & turn-key
- Integrated dispersion pre-compensation

APPLICATIONS

- Multiphoton microscopy
- Neuroscience
- Photopolymerization
- Ophthalmology
- OPO pumping

BIOLIT 2

The **Biolit 2** is a compact, air-cooled femtosecond laser designed for multiphoton microscopy, biophotonics and other non-linear optics applications.

The industrial-grade device is exceptionally robust, maintenance-free and affordable. A combination of ultra short (typ 60 fs) and clean pulses, superior beam quality and lower (compared to solid-state siblings) repetition rate enables unparalleled multiphoton imaging while preserving the object.

SPECIFICATIONS

| Model | Biolit 2 | Biolit 2 SH ¹⁾ |
|---|--|---------------------------|
| Central wavelength | 1050 ± 5 nm | 525 ± 5 nm |
| Average power | > 2 W | > 400 mW |
| Pulse duration | < 70 fs (60 fs typ.) | |
| Pulse duration stability ²⁾ | < ± 5 fs | |
| Pulse strehl ratio | > 0.9 | |
| Tunable dispersion pre-compensation ³⁾ | -10 000 fs ² ... +500 fs | N/A |
| Pulse repetition rate ⁴⁾ | 15, 20, 30 or 40 MHz | 15 or 20 MHz |
| Tuneable pulse repetition rate (optional) ⁵⁾ | 1 – 40 MHz | |
| Analog power control | 1 – 100 % | |
| Analog power control bandwidth | > 10 kHz – standard >300 kHz – optional | |
| Peak power | > 1 MW @ 20 MHz | >200 kW @ 20 MHz |
| Beam quality | $M^2 < 1.2$ | |
| Beam circularity ⁶⁾ | > 0.9 | |

| Model | Biolit 2 | Biolit 2 SH ¹⁾ |
|--|--------------------------------|---------------------------|
| Beam diameter (1/e ² level) | 1.5 ± 0.3 mm | 1.2 ± 0.3 mm |
| Beam pointing (pk-to-pk) ²⁾ | < 50 µrad | |
| Beam pointing vs temperature | < 25 µrad/°C | |
| Pulse energy stability (RMS) | < 1 % | < 2 % |
| Power stability (RMS) ²⁾ | < 1 % | < 2 % |
| Warm up time (cold start) | < 10 min | |
| Available control interfaces | USB, CAN | |
| Powering requirements | 100 ... 240 V AC, 47 ... 63 Hz | |
| Operating temperature | 15 – 35 °C | |
| Humidity | non condensing | |
| Transportation/storage temperature | - 20 – +70 °C | |
| Dimensions: | | |
| Laser head (L × W × H) | 334 × 217 × 139 mm | 334 × 211 × 165 mm |
| Control unit (L × W × H) | 449 × 370 × 140 mm | 449 × 370 × 140 mm |
| Umbilical length | 3 ± 0.3 m | |
| Colling: | | |
| Laser head | air (passive) | |
| Control unit | forced air (fans) | |

¹⁾ Biolit-2 SH model also has infrared (1050 nm) output with the same specifications as standard Biolit laser. Both outputs are available simultaneously.

²⁾ Measured during 24 h operation after 10 min warm-up.


³⁾ Equivalent of 100 mm of SF10 glass. Even higher dispersion (up to 30'000 fs²) of the external optical system can be pre-compensated in the factory on request.


⁴⁾ Factory preset. Other repetition rates are available on request. Please inquire for more details.

⁵⁾ Output repetition rate can be described by formula $RR = RR_0/N$, where RR_0 is fundamental repetition rate and N is integer number. Output power is dependent on both RR and RR_0 . For power dependence on the repetition rate please contact LITILT.


⁶⁾ Defined as the worst case ellipticity along the z-scan ($\pm 5 \times L_{\text{Rayleigh}}$) of the beam.

World patented technology: US10038297, JP6276471, EP3178137, CN106575849.







Output spectrum from **Biolit 2** laser




Measured autocorrelation function of the pulses from **Biolit 2** laser




Long term power stability of **Biolit 2** laser (at 1050 nm)



Beam diameter dependence on propagation distance (z-scan) of **Biolit 2** laser and M^2 fit




Drawing of **Biolit 2** laser head (in mm)



Drawing of **Biolit 2 SH** laser head (in mm)





Drawing of **Biolit 2** power/control supply (in mm)